



- Reduced size and weight (13 grams)
- TempVerified updated software (Alarm: 39.7c plus 6 hours)
- Automatic alarm reset
- Treatment success validation
- Designated installation location
- New patent pending bio-degradable plastic compound
- Flexible tag feature
- Probe retention rotational feature
- Weather and water proof
- Clinically accurate (undisturbed) temperature measurement

## FEVERTAGS & E-TAG

REDUCED
SIZE AND WEIGHT

UPDATED SOFTWARE

SPECIFIC INSTALLATION LOCATION

WATER PROOF ENCAPSULATION

**EASY TO INSTALL** 

