



HEATED THEARAPY

- TAIWAN AND CHINA PATENTED
- Carbon Fiber Heating Technology
- Infrared Heat for Pain Relief & Other Health Benefits









STAYWARM®

Movable Rechargeable Cordless Knee / Back / Shoulder / Elbow / Wrist / Ankle Heat Wrap from SUNGROUP

- **Relax Muscle Spasms**
- **Reduce Chronic Pain**
- **Improve Circulation**

System

- Generate Soothing Warm
- **Boost Metabolism & Immune**

Fast Heat Up

No Electro Magnetic Field

FAR Infrared

Save Money

Hand Wash

3 Level Heating

Taiwan MEDICAL DEVICE LICENSES APPROVED

Taiwan and China Medical Licences Approved TFDA Produced under norms and conditions for International







MELTS AWAY MUSCLE TENSION AND FATIGUE.

Our knees and Back take on a great deal of trauma which can lead to serious pain and discomfort. Our Heated Knee and Back Wrap helps to soothe the pain and increase blood flow without the use of potentially harmful medicines.

If you are ready for relief, you are ready to experience the benefits that our Heated wrap has to offer. Our system delivers therapeutic and targeted heat directly to your muscles and tissues where you need it most. Regain your mobility and active lifestyle again with the At-Home Heated Wrap



PREVENTION & Management From Sports Injuries

Playing sports is good for your health, but sometimes you can injure yourself. Prevention of training injuries is far easier than fixing them afterwards.

There are so many reasons sports injuries can be caused by, for example, not warming up properly before exercising and pushing yourself too hard or using inadequate equipment. According to the American College of Sports Medicine, up to 50 % of all athletic injuries can be avoided if athletes take appropriate preventive action.

SOLUTION FOR SYMPTOMS THAT CAUSE PAIN AND DISCOMFORT AND PRE INJURY MANGEMENT.

Timely use of Stay warm Wrap after any intence activity or sports can prevent muscle and Ligaments from fatigue, Long time negligence taken could cause serious bone or muscles injuries.

Our FIR Heat Therapy Wrap relaxes your muscles, tendons, and ligaments with soothing heat to increase your flexibility and blood circulation thereby relieving the sensation of pain.


How does it work? Staywarm Heat Therapy Wraps are constructed with our patented pending carbon fibre heating system that uses tiny un noticeable conductive fibers that provide instant heat.

The LED Color indicator on the product clearly indicates the required setting, which makes it easy to identify the power setting on the device. With similar products the battery pocket is stationary and in some cases very inconvenient, but Staywarm Therapy has incorporated our new flex pouch battery pocket that allows you to place the battery anywhere on the wrap. This allows each individual user to place the battery in the locations that best suits them and their use.

The Staywarm Heat Therapy Wraps heat through conduction. This means that the heat flow is transferred into your joints and muscles due to direct contact.

By using conduction, we can effectively heat targeted areas using a rechargeable 7V Lithium battery.

Staywarm Heat Therapy Wraps provide hours of soothing warmth to help rejuvenate sore joints and tired muscles.



DEEP PENETRATION THERAPY

FAR INFRARED RAY TECHNOLOGY

F.I.R. Technology transforms light energy emitted by the sun into therapeutic heat energy that stimulates the body's natural healing capabilities and lymphatic movements to get rid of toxins and oxidative stress.

HEAT THERAPY

Heat therapy is the most recommended non-invasive, non-pharmaceutical form of pain relief and rehabilitation by physical therapists and medical professionals alike.

CONVENIENCE

The revolutionary design brings the power of Far Infrared lamps seen at the doctor's office into the convenience of your home!

SAFE LOW-VOLTAGE

Our heat therapy products are low-voltage and produce virtually no electromagnetic fields (EMF) for your safety.

Taiwan MEDICAL DEVICE LICENSES APPROVED

SPECIFICATION

Safe, Far Infrared Ray (FIR) Heat Automatic 5-Minute Safety Timer Personalized Temperature Controller

3D TAILORED. ADJUSTABLE. FORM-FITTING.

3D tailored and anatomically designed garment positioned

directly on the elbow back area to fully relax tense muscles, ligaments, and tendons.

Rechargerble battery

We use the best Lithium ion rechargeble Cell from Samsung and LG.



Please use carefully, may cause serious burns, consyour physician if pain persists.



The unattended use of this heated therapy wrap by children or incapacitated person may be dangerous.

Innovation for Your Well-Being









