





# Sterilization Centers



# • E155 Sterilization Center

Luxurious style. Quality performance.



# Hands-free Operations

Color-coded pneumatic switches to operate hands-free waste drawer (green), unsterile upper storage doors (red), and hands-free tower doors (blue).

#### Sterile

Red light storage for soiled instrument cassettes and sanitized instrument storage with blue light to indicate items have been processed.

### Safe Environment

Pneumatic hands-free waste drawer with sharps easily manages waste materials based on safety guidelines. The push-out opening mechanism allows the user to open and close the drawer without any cross-contamination.

## Accessibility

Tapmaster<sup>™</sup> foot control can be found underneath the sink for hands-free operations.

















# E120 Sterilization Centers

Embodies a customizable form and premium function.



#### Sterile

Pull-out rotating soft-closing instrument tower with blue light.

### Prep

Unsterile upper instrument cabinet with glass shelves and red light.

# Hands-free Operations

Tapmaster<sup>™</sup> foot control can be found underneath the sink for hands-free operations.

# Convenience

Pull-out waste container and full-extension soft close tandembox drawer.

# Upgrades

Options include task lights, sharps disposal container, sink upgrades.

















# D96 Sterilization Centers

Incorporates a modular design with customizing components to deliver simplistic function.



### Sink Usability

Tapmaster<sup>™</sup> kick pedal can be found underneath the sink for hands-free operations

# Workspace Components

Standard laminate countertop with drop-in stainless steel sink and chrome faucet with soap dispenser.

# Steri-Tower

Sterilization tower utilizes wrap drawers with slide-out shelves to keep tools accessible and within reach.

# Clean Atmosphere

Pull-out trash module keeps waste out of sight for a cleaner look.

Belmont Equipment 101 Belmont Drive Somerset, NJ 08873 (800) 223-1192 www.belmontequip.com Koken Manufacturing Company, Inc. A Division of Takara Belmont USA, Inc. 1631 Martin Luther King Drive St. Louis, MO 63106 T: (800) 325-7373 F: (314) 446-7088

BelmontEquipUS







9

in



Takara Company Canada Ltd. 2076 South Sheridan Way Mississauga, ON L5J 2M4, Canada (800) 268-5351 dental.takarabelmont.ca

