



Alfonso Piñeyro, DDS

Dr. Piñeyro received his dental degree from the Universidad Autonoma de Guadalajara. He then completed a two year advanced education in general dentistry before completing his specialist training in prosthodontics at the University of Rochester Eastman Dental Center.

Dr. Piñeyro is currently an affiliate instructor at the Department of Restorative Dentistry at the University of Washington. He is the past president of the esteemed Washington State Society of Prosthodontists and is also a member of the American College of Prosthodontists and the American Prosthodontics Society.

Dr. Piñeyro has authored several publications in peer reviewed dental journals, chapters in textbooks and is currently involved in research focusing on dental implants. Dr. Piñeyro practices full time in a practice limited to fixed, removable and implant prosthodontics.

He has lectured both nationally and internationally and has presented over 250 lectures for multiple dental study clubs and dental associations. He has been invited to speak to large national dental symposiums which include: Seattle Study Club Symposium, American Academy of Periodontology, American Academy of Implant Dentistry, American Academy of Fixed Prosthodontics, Academy of Osseointegration, American College of Prosthodontists, Academy of LDS Dentists and The 12th International Symposium on Periodontics and Restorative Dentistry.

The Pineyro Arch Kit

The latest innovation in implant maintenance from the creators of the Montana Jack®

For more information visit www.pdt dental.com/pineyro



Most **ergonomic** handle with knurling to instrument shank providing comfort with less hand fatigue and pinch.

With a solid resin handle weighing only **13 grams**, PDT has the **lightest**, most **tactile sensitive** instruments on the market.

Anatomically color-coded for **fast** set-ups and **efficient** pick-ups.

Instruments retain their **sharpness** due to the hardness of the steel.

VISIT PDT.DENTAL/FULLCATALOG

to view our full line of instruments including: curettes & scalers, sharpening kits, bone grafting instruments, composite & restorative instruments, crown removing pliers, elevators, extracting forceps, implant impression trays, needle holders, periostomes,scalpel handles, scissors, tweezers, and more!

PDT, Inc. - Paradise Dental Technologies
8275 Highway 10 W ♦ Missoula, MT 59808
Phone: 800-240-9895 ♦ Fax: 406-626-4550
customer.care@pdt dental.com ♦ pdt dental.com

The Pineyro Arch™ Kit

Patent Pending


The first set of instruments specifically designed for the implant debridement of full-fixed implant restorations.



Allows for the treatment of full-fixed hybrid patients without needing to remove the prosthesis.



PARADISE DENTAL TECHNOLOGIES




Pineyro Arch Ti 1

Anterior Arch - to debride the anterior placed implants for a Full-Fixed Hybrid, especially lingual area of lower anterior teeth, where calculus builds up most frequently




Pineyro Arch Ti 2

Posterior Arch - to debride the posterior placed implants for a Full-Fixed Hybrid, the design of the shank allows for proper access



Pineyro Arch Ti 3

Mesial/ Distal Arch - to debride the posterior hard to reach mesial/ distal areas on the implants for a Full-Fixed Hybrid



Pineyro Arch Ti 4

Specialty Arch - or all implants prior to taking impressions -clean off abutments for single to hybrid cases prior to cementation

No more removing screws for hygiene appointments
 Save Time ♦ Save Frustrations ♦ Save Calls To Your Implant Rep



www.pdt dental.com/pineyro

R930 - Pineyro Full Arch Kit
 Includes T005A-F Fliptop Cassette

Made of solid titanium that won't scratch or leave debris

Used in up and down or side to side motions

