

Callisto™

Portable audiometry & REM



Redefining portability

- Audiometry, REM, and Visible Speech Mapping in one box



Interacoustics®

leading diagnostic solutions

Callisto™

Portable audiometry & REM

Fitting into a modern practice



Visible Speech Mapping (optional)

- Intuitive – combines advanced testing and counseling
- Conforms to other Interacoustics Suite Software reducing the learning curve
- Mix and match Audiometry, REM and VSP
- Tailored customizable test routines
- Merge Video Otoscopy into the suite
- NOAH/OtoAccess™ compatible



Bring the clinic to your client

- Ideal for travelling and in-home visits
- Small footprint and lightweight (1.25 lbs)
- Powered via the PC USB
- Optional dedicated carrying bag on wheels
- Leave Callisto in the bag while testing



Carrying bag

Carrying bag inside



"Inspired by advanced Audiometry and Real-ear systems, but wrapped in a small package"

Audiometry


- Customizable test setups, printouts and report templates
- Visual and auditory counseling tools (Speech Banana, Hearing Loss Simulator, Master Hearing Aid)
- Digitized word lists from your PC hard-drive
- Audiogram historical comparison overlay
- High Frequency capability and Digitized QuickSIN™ (optional)



Audiometry

REM

- Visible Speech Mapping for combined verification and counseling (optional)
- Multiple speech stimuli in several languages and live voice
- True 'On-Top' mode for fast interactive HA fittings
- New NAL-NL2 and DSLv5
- Open Fit capable
- Directionality testing
- See data in graph or table mode
- Quickly switch between Response or Gain views
- In-the-ear monitor
- RECD capability



Real Ear Measurement (REM)



Interacoustics®

leading diagnostic solutions

Interacoustics – the wise choice

With over 40 years of experience, Interacoustics is dedicated to supplying its customers with the best possible solutions for their audiologic needs. This is accomplished by maintaining a continuous dialogue with healthcare professionals working in all sectors of audiology. Our equipment meets the highest possible engineering standards and we provide design expertise that can only come from close contact with clinical practice.

Solutions on every scale

Designing equipment for every size of clinic in so many countries puts us in the unique position of being able to offer solutions that fit your requirements exactly. Audiometry, tympanometry, electrophysiology, hearing aid testing, balance investigation are all within our scope and can be integrated to suit your needs.

Designed for diagnostics efficiency

We design equipment to make testing and interpretation easier. This means better interfaces, well designed screen layouts, printed reports and interaction over networks with databases and electronic records systems. In most cases, you can configure the settings and layout yourself.

Support worldwide

The Interacoustics name is not only your guarantee of quality and functionality, but also for support. We operate in over 100 countries worldwide through a well coordinated network of distributors and service centres to ensure that you receive total support and service.

Safety standards:	IEC 60601-1, Applied parts Type B, UL 60601-1, CAN/CSA-C22.2 No. 601.1-M90, IEC 60601-1-1
Computer communication:	USB interface, compatible with USB1.1 or later
PC minimum requirements:	General: 2 GHz Intel Core 2 Duo CPU (or better) with at least 1GB ram (preferably 2GB), DVD drive, USB connector. Display: minimum resolution of 1024x768 with hardware accelerated DirectX/Direct3D graphics card. Disk Space: At least 1.5 GB available space
Operating system:	Windows 7 (32 and 64 bit), Windows VISTA (32 bit) and Windows XP (SP2 or later and compatible)
Database support:	OtoAccess™ and NOAH 3.7.x compatible office systems or later releases
Construction:	Plastic cabinet
Power:	USB-powered with internal "power boost" rechargeable battery and load balancer
Weight:	565 g / 1.2 lbs.
Dimensions:	212 x 121 x 44 mm / 8.3 x 4.8 x 1.7 inches

Sales and service in your area:

Interacoustics USA

Phone: 1-800 947 6334 · Fax: 1-952-903-4200
E-mail: info@interacoustics-us.com
7625 Golden Triangle Drive, Eden Prairie, MN 55344
Web: www.interacoustics-us.com

Products in this group:

- Callisto™ hardware
- REM440 Real Ear Measurement module
- AC440 Audiometry module

Optional special tests for AC440:

- High Frequency audiometry (HF440)
- Speech from Hard-drive (SFH440),
- SISI test
- Master Hearing Aid (MHA440),
- Hearing Loss Simulator (HLS440)
- QuickSIN™

Optional special tests for REM440:

- VSP440 Visible Speech Mapping and Counselling module

Related products:

- Affinity^{2.0} Hearing Aid Analyzer
- Equinox^{2.0} Clinical PC-based audiometer



Add an overlay image that matches your clinic

83020801-US -1 - 09/2011

Read more here:
www.interacoustics-us.com/Callisto



Interacoustics®

leading diagnostic solutions