

PRT-TEST

DRY EYE TESTING IN JUST 15 SECONDS

PHENOL RED THREAD TEST

The “Colour-Change” Dry Eye Test

- Can Provide a Realistic Measure of the Resting Tear Volume.
- **Less Reflex Tearing** compared to other tests
- **Painless & Comfortable** for the end-user
- Recent studies showing a **Good Correlation** with the **Schirmer’s Test**
- Can **Combine with the Schirmer’s Test** to give more effective & accurate results

Optimum Design & Premium Quality

- 100 Sterile Threads wrapped in Transparent pouches
- Produced using the Highest Grade Material
- Designed in the UK & Manufactured under GMP

AS RECOMMENDED BY



tfosDEWS II



For more Information Contact:

Entod Research Cell UK Ltd.

15 Tottenham Lane, Hornsey, London N8 9DJ, UK.

Tel: (+44) 0208 144 0963

Email: contact@ophthalmicsuk.com