

AC40

Clinical Audiometer

RELEASE
FOR SALE 2012



The New AC40 True Hybrid Audiometer - standalone and PC-audiometer in one solution!

- **NEW** Pediatric noise stimuli
- **NEW** New Masking Help and Auto Masking features (Diagnostic Suite 2.0)
- **NEW** TEN (HL) test for diagnosing dead regions in the cochlea
- **NEW** Patient counseling tools
- **NEW** Direct HDMI output for external monitors

The new AC40 clinical hybrid audiometer includes all the advanced features and functions needed by a modern high-end clinic. It provides the power and capabilities of a PC-based instrument, without always being dependent on a computer. The AC40 concept incorporates clinical efficiency & ease of use without compromise.

More highlights

- True two channel with full clinical test battery
- Flexible reporting tools – combine results with the Titan middle ear analyzer
- Large 8.4" articulating color display
- Exceptional power amplifier
- Fully EMR compatible (XML and PDF)
- Store results in Noah or OtoAccess™
- Internal storage for 1,000 patients / 50,000 test sessions
- Integrated pre-recorded speech
- Direct print options



Interacoustics®


leading diagnostic solutions

AC40

Clinical Audiometer

Define
multiple user
test protocols
on the fly

Built-in
wave files




Speech table screen. Also available as speech graph



Combine test results with the handheld impedance platform. Titan is available as a diagnostic or clinical tympanometer including optional automated ABR & OAE modules.

New! Counseling overlays:

- Phoneme examples
- Sound examples
- Speech banana
- Severity indications



New software tools provide an added dimension in patient counseling

The Interacoustics True Hybrid Concept

THE TRUE HYBRID CONCEPT ENABLES YOU TO:

1. Operate the AC40 as stand alone
2. Operate the AC40 via computer
3. Operate the computer via the AC40

Both the AC40 clinical audiometer & the AD629 diagnostic audiometer as well as the Titan diagnostic/clinical middle ear analyzer are true hybrids

true
hybrid

Sales and service in your area:

Check www.interacoustics-us.com
for updates....

Interacoustics USA

Phone: 1-800 947 6334 · Fax: 1-952-903-4200
E-mail: info@interacoustics-us.com
7625 Golden Triangle Drive, Eden Prairie, MN 55344
Web: www.interacoustics-us.com



Interacoustics®

leading diagnostic solutions